



Compact Kit Technicians' Choice

Zirconia Contouring Kit — For Technician

Special compact kit for dental laboratory.

It is mainly used for fast contouring and polishing of zirconia ceramic materials. It is also suitable for difficult-to-polish parts such as the occlusal surface, which is more precise and convenient. It can achieve high gloss without glazing.



				
Order No.	TB11C	TG03	TG10	HT06C
Size(mm)	22*2	4.5*13	13*2	20*2
Rotation speed(rpm)	8000-15000	8000-15000	8000-15000	7000-12000
Material	All-ceramic grinding head coarse grinding wheel shape	All-ceramic grinding head fine sand conical shape	All-ceramic grinding head fine sand flat head	Whirlwind wheel finisher
Usage	Surface shaping	Surface shaping	Surface shaping	Marginal polishing

				
Order No.	FD2222B	MP18	MD15	FK30
Size(mm)	22*0.22	1.8*45	1.5*45	3*45
Rotation speed(rpm)	7000-12000	8000-15000	8000-15000	8000-15000
Material	Elastic emery sheet	Multi-coated emery burs	Multi-coated emery burs	Multi-coated emery burs
Usage	Crown and bridge gap dealing and shape modification	The groove adjusting and interdental trimming	The groove adjusting and interdental trimming	The groove adjusting and interdental trimming



Compact Collocation Best Chairside Assistant

Zirconia Polishing Kit—— For Clinical Use

It is mainly used for fast contouring and polishing of zirconia ceramic materials, and is suitable for difficult-to-polish parts such as occlusal surface. Simple and efficient operation. Only polishing is required to complete the restoration, which is suitable for clinical operation.



						
Order No.	HT10M	HT03M	HG17M	HG03M	HG17F	HG03F
Size(mm)	12*2	4*10	17*1.6	4*13	17*1.6	4*13
Rotation speed(rpm)	8000-12000	8000-12000	7000-15000	7000-15000	7000-15000	7000-15000
Material	Wheel refiner	Cone refiner	Diamond polishing grinding head cyclone wheel shape	Diamond polishing grinding head bullet type	Diamond polishing grinding head cyclone wheel shape	Diamond polishing grinding head bullet type
Usage	Marginal polishing	Marginal polishing	Pre-polishing	Pre-polishing	high-gloss polishing	high-gloss polishing

Note: The rotation speed of the polishing tools should be operated strictly according to the instruction.

The grinding head of coarse sand is sharp, and the grinding is fast but the surface is not smooth as fine sand.

The cyclonic wheel shaped product cannot be polished with too much pressure, and it needs to be polished in a clockwise direction, otherwise it is easy to break.