

Bloomden Universal Nano Hybrid Composite Resin

Fast curing

High wear resistance

Low polymerization shrinkage


• Better trust Better smile



Bloomden Universal Nano Hybrid Composite Resin


Bloomden Universal Nano Hybrid Composite Resin

Nanofiller, used for anterior and posterior restorations. Low polymerization shrinkage, high strength, and good operation. Available in 5 colors, easy to use.



Application

- Direct restoration of various types of holes in anterior and posterior teeth
- Used in combination with glass ionomer for "sandwich" fillings
- Cusp reshaping
- Core Build Up
- Fixed splints
- Indirect restoration of anterior and posterior teeth: veneers



Features

- Low polymerization shrinkage, greatly reducing the pulling force and micro-tension on tooth tissue, effectively reducing postoperative sensitivity.
- Good wear resistance, excellent anti-fracture and strength, making the restoration more durable
- Good operation, easy to fill and easy to make the shape before curing, saving working time
- Fast curing in 10-20 seconds, with a curing depth of up to 2.5mm

Shade

A1,A2,A3,A3.5,B1

Specification

4g/stick