



Bloomden Opacity Liquid

Predictable Aesthetics
through Controlled Opacification



Bloomden Opacity Liquid

Specially formulated for zirconia, it delivers effective masking while preserving material translucency, achieving the ideal balance between opacity and vitality.



Reliable Masking & Uniform Base

Effectively masks dark base including tetracycline-stained teeth, fluorosis, and metal post/cores. Creates a uniform, light base layer that prevents show-through, ensuring true-color restoration.



True Color Rendition & Controllable Value

Ensures accurate color reproduction for true "what-you-see-is-what-you-get" results. Supports high-value restorations unrestricted by tooth substrate color.



Effortless Application & Precision Control

Low-viscosity liquid enables thin, even application with minimal brush effort. A single thin layer provides full coverage without compromising seating.

