

Bloomden Diamond Paste

- Effective
- Easy polishing
- Easy Clean up

• Better trust Better smile



Bloomden Diamond Paste

Advantage

Bloomden Diamond Paste is a solid polish in a jar that can be used to polish high gloss polished ceramic materials such as lithium disilicate, zirconium oxide, veneers and hybrid ceramics, as well as high performance polymers and veneer composites.



Effective & Ease-of-use

Bloomden Diamond Paste is engineered for optimal performance with a user-friendly application. Its formulation guarantees effective results with minimal effort, making it the preferred choice for dental professionals who value efficiency and high-quality outcomes.



Easy polishing in occlusal anatomy

This paste excels in polishing the detailed surfaces of occlusal anatomy. The fine abrasive particles are specifically designed to allow precise and thorough polishing in the most intricate areas, ensuring a smooth and natural finish for every part of the restoration.



Skip the Glaze without the seashell effect

Most Diamond polish leaves an opalescent / seashell effect. Bloomden Diamond Paste disrupts the convention of glazing all your crowns. With it, you can mill & polish your crowns to a life like finish without the need for the final glaze.



Dust-free & Easy Clean up

Made of natural materials, Bloomden Diamond Paste delivers the perfect combination of diamond abrasive media delivered in a premium, dust-free, water soluble binder that is easy to clean.

Component name	Ingredient	Content%	CAS
Diamond powder	C	1-8%	7782-440-3
Aluminium oxide	Al ₂ O ₃	1-8%	1344-28-1
glycerin	C ₃ H ₈ O ₃	≥70%	56-81-5